

Optional Carrying Case and Portable Stand for HandyRef-K/HandyRef

1 Carrying case for HandyRef series with portable stand

- Holds both HandyRef-K/HandyRef and portable stand
- Portable stand is included with the case

Dimensions: 505 (W) x 265 (D) x 390 (H) mm
Mass: 8.5 kg (including portable stand)*

The case holds	
HandyRef-K/HandyRef	✓
Portable stand	✓

HandyRef-K/HandyRef sold separately



2 Carrying case for HandyRef series

- Holds HandyRef-K/HandyRef only
- Does not hold portable stand

Dimensions: 444 (W) x 197 (D) x 324 (H) mm
Mass: 3.7 kg*

The case holds	
HandyRef-K/HandyRef	✓
Portable stand	

HandyRef-K/HandyRef sold separately



3 Carrying case for portable stand

- Dedicated for portable stand alone
- Portable stand is included with the case
- Does not hold HandyRef-K/HandyRef

Dimensions: 497 (W) x 120 (D) x 411 (H) mm
Mass: 5.5 kg (including portable stand)*

The case holds	
HandyRef-K/HandyRef	
Portable stand	✓

HandyRef-K/HandyRef sold separately



*Total weight of the case and included items.

Specifications may vary depending on circumstances in each country.
Specifications and design are subject to change without notice.

