



Injector for Nex-Acri™
Nex-IJ




THE ART OF EYE CARE

Handpiece

Product name	Nex-IJ
Model	IJ-2

Do not use this device to implant intraocular lenses other than the NIDEK Nex-Acri AA 1P, Nex-Acri AA or Nex-Acri intraocular lens. Any use other than those specified may cause unexpected malfunction or adverse device effects and complications.



Cartridge

Product name	Nex-IJ		
Model	Type 4C	Type 3C	Type 2C
Lens	Nex-Acri AA 1P (NS-60YG) Nex-Acri 1P (NS-60G)	Nex-Acri AA (N4-18YG)	Nex-Acri AA (N4-18YG) Nex-Acri (N4-18B)*
Incision Size	2.4 mm	2.4 mm	2.8 mm

*Applicable lens diopter range is 1.0 to 28.0 D.

This device is designed only for use in combination with the NIDEK Nex-IJ handpiece. Do not use this device alone or use it with other medical devices.



Type 4C



Type 3C



Type 2C



CAUTION

To use this product safely, please read instructions carefully before use.

Brochure and listed features of the device are intended for non-US practitioners. Specifications may vary depending on circumstances in each country. Specifications and design are subject to change without notice.



HEAD OFFICE
(International Div.)
34-14 Maehama,
Hiroishi-cho, Gamagori,
Aichi 443-0038, JAPAN
TEL: +81-533-67-8895
URL: www.nidek.com
[Manufacturer]

TOKYO OFFICE
(International Div.)
3F Sumitomo Fudosan Hongo
Bldg., 3-22-5 Hongo, Bunkyo-ku,
Tokyo 113-0033, JAPAN
TEL: +81-3-5844-2641
URL: www.nidek.com

NIDEK INC.
2040 Corporate Court,
San Jose, CA 95131, U.S.A.
TEL: +1-408-468-6400
+1-800-223-9044
(US Only)
URL: usa.nidek.com

NIDEK S.A.
Ecoparc,
9 rue Benjamin Franklin,
94370 Sucy En Brie,
FRANCE
TEL: +33-1-49 80 97 97
URL: www.nidek.fr

NIDEK TECHNOLOGIES S.R.L.
Via dell'Artigianato,
6/A, 35020 Albignasego (Padova),
ITALY
TEL: +39 049 8629200/8626399
URL: www.nidektechnologies.it

NIDEK (SHANGHAI) CO., LTD.
Rm3205, Shanghai Multi
Media Park, No.1027 Chang
Ning Rd, Chang Ning District,
Shanghai, CHINA 200050
TEL: +86 021-5212-7942
URL: www.nidek-china.cn

NIDEK SINGAPORE PTE. LTD.
51 Changi Business Park
Central 2, #06-14,
The Signature 486066,
SINGAPORE
TEL: +65 6588 0389
URL: www.nidek.sg